



# INFUSION SETS

For most pumps, the insulin in the device needs help getting to your body. That's where an infusion set comes in. The flexible tubing connects to a pump reservoir and carries the insulin from the pump to the needle or cannula beneath your skin.

COMPANY	BRAND	PUMP COMPATIBILITY	NEEDLE/CANNULA	TUBING*	INSERTION METHOD	
<b>ANIMAS CORP.</b>	Comfort and Comfort Short	Compatible with all Luer-lock pumps	20° to 45° insertion angle; 25-gauge plastic cannula; 13-mm (Comfort Short) and 17-mm (Comfort) cannula lengths	23, 31, and 43 in.; at site	Manual	
	Contact Detach		90° insertion angle; 29-gauge stainless steel needle; 6- and 8-mm needle lengths	23 and 43 in.; away from site		
	Inset		90° insertion angle; 25-gauge plastic cannula; 6- and 9-mm cannula lengths	23 and 43 in.; at site		Via built-in disposable insertion device
	Inset 30		30° insertion angle; 25-gauge plastic cannula; 13-mm cannula length			
<b>MEDTRONIC DIABETES</b>	MiniMed Mio	Compatible with Medtronic pumps only	90° insertion angle; 25-gauge plastic cannula; 6- and 9-mm cannula lengths	18, 23, and 32 in.; at site	Via built-in disposable insertion device	
	MiniMed Quick-Set	Compatible with Medtronic pumps and pumps that have a Luer-lock connector		18, 23, 32, and 43 in.; at site	Manual or via MiniMed Quick-serter insertion device	
	MiniMed Silhouette			20° to 40° insertion angle; 25-gauge plastic cannula; 13- and 17-mm cannula lengths	Manual or via MiniMed Sil-serter insertion device	
	MiniMed Sure-T			90° insertion angle; 29-gauge stainless steel needle; 6-, 8-, and 10-mm needle lengths	18, 23, and 32 in.; away from site	Manual

COMPANY	BRAND	PUMP COMPATIBILITY	NEEDLE/CANNULA	TUBING*	INSERTION METHOD
SMITHS MEDICAL	Cleo 90	Compatible with all Luer-lock pumps	90° insertion angle; 25-gauge plastic cannula; 6- and 9-mm cannula lengths	24, 31, and 43 in.; at site	Via built-in disposable insertion device
SOOIL DEVELOPMENT	Easy Release	Compatible with Sooil pumps only	90° insertion angle; 27-gauge stainless steel needle; 7- and 9-mm needle lengths	27 and 43 in.; away from site	Manual
	Soft-Release O		90° insertion angle; 26-gauge plastic cannula; 6- and 9-mm cannula lengths	24, 31, and 40 in.; at site	
	Soft-Release ST		15° to 30° insertion angle; 24-gauge Teflon cannula; 19-mm cannula length	27 and 43 in.; at site	
TANDEM DIABETES CARE	AutoSoft 30	Compatible with Tandem pumps only	30° insertion angle; 25-gauge plastic cannula; 13-mm cannula length	23 and 43 in.; at site	Via built-in disposable insertion device
	AutoSoft 90		90° insertion angle; 25-gauge plastic cannula; 6- and 9-mm cannula lengths		
	TruSteel		90° insertion angle; 29-gauge stainless steel needle; 6- and 8-mm needle lengths	23 and 32 in.; away from site	Manual
	VariSoft		20° to 45° insertion angle; 25-gauge plastic cannula; 13- and 17-mm cannula lengths		

**\*At site:** Tubing disconnects at the infusion site, leaving only the needle or cannula and adhesive in place. **Away from site:** Tubing disconnects a few inches away from the infusion site; the needle or cannula, adhesive, and a short tail of tubing remain.